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In Situ Monitoring of Curing Reactions with Mid-Infrared **Laser Spectroscopy**

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Fast spectroscopic monitoring is a proven tool in R&D and production environments, where it is used to gain insight to chemical reactions. This application note showcases how the IRis-F1 dual-comb spectrometer can be used to monitor curing reactions, where it can provide invaluable insight on the performance of coatings and adhesives.

Infrared spectroscopy is an ideal measurement technique

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for this application in many ways as it is capable of performing non-invasive in-situ measurements that provide direct information about the structure and bonding of the sample. However, many high-performance adhesives possess fast curing rates and the speed of off-the-shelf spectroscopic techniques is the limiting factor in their analysis.

The IRis-F1 is based on quantum cascade laser (QCL) frequency comb technology, which allows for measurements on milli to microsecond timescales to be carried out.

Experimental

The experimental setup is depicted in Figure 1. Measurements were conducted with a single-reflection diamond Golden Gate ATR accessory (Specac). The sample is a UV-triggered adhesive, which was kindly provided by DELO Industrial Adhesives (Germany). It possesses an infrared band at $1,613\text{cm}^{-1}$ when uncured and a band at $1,638\text{cm}^{-1}$ after curing. For each measurement a background file was initially acquired after cleaning the ATR crystal, followed by a rapid series of measurements spanning before, during, and after curing.

The QCL used provides a spectral output of $1,585 - 1,645\text{cm}^{-1}$. The spectrometer was set to a resolution of 4cm^{-1} and the acquisition time was 1ms , with a repetition rate of 40Hz . At these settings, continuous acquisitions can be acquired for ca. 30 seconds before the buffer fills and repetition rate is reduced to the order of 1Hz . Illumination was provided using a Hamamatsu LC8 mercury lamp, via an optical fibre (Figure 1, right). After the measurement is started on the spectrometer, the lamp is turned on to initiate the curing reaction and the resulting spectra are recorded as a time-series.

Results and discussion

Figure 2 shows the raw spectra acquired in time-increments of 25ms . A good signal to noise ratio is obtained even at short integration times. It is clear that the band at ca. $1,610\text{cm}^{-1}$ decays rapidly, while a new band grows at ca. $1,640\text{cm}^{-1}$. The fringes observed in the raw spectra are expected when using


图1: (左) 带有Golden Gate ATR附件的IRis-F1光谱仪; (右) 安装在钻石反射元件上方的UV光源。
Figure 1: (left) The IRis-F1 spectrometer with the Golden Gate ATR accessory; (right) the UV light source mounted above the diamond reflection element.


图2：以25ms的时间分辨率测量的固化过程的光谱；使用FT滤镜去除相同的光谱，以去除ATR棱镜上的干涉条纹。
Figure 2: Spectra of the curing process measured with a time-resolution of 25 ms; the same spectra with an FT filter to remove interference fringes from the ATR prism.

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a laser-based instrument with an ATR prism. They are caused by a change in the refractive index of the sample (i.e. glue vs air) changing the penetration depth in ATR spectroscopy. Fringes are typically not observed with traditional light sources, as there the coherence length is short enough to be unaffected by fringes on this length-scale. Thanks to their well-defined periodicity, they can be removed with a simple Fourier transform (FT) filter (Figure 2).

Difference spectra can be generated by subtracting a spectrum of the uncured glue (defined as a negative-time spectrum) from each of the subsequent positive-time spectra. Figure 3 (left) shows the difference spectra, which highlight the changes produced by the curing reaction. Negative bands indicate that a feature has disappeared, while positive bands indicate that a new feature has grown in. Here, it is evident that not one, but two new bands has emerged, at ca. 1,638 and 1,624 cm^{-1} . It should be noted that no filtering was applied to the spectrum in Figure 3 (left).

The kinetic traces at two points in the spectra clearly show the time-point where the lamp was turned on. By fitting the growth and decay of these points, the kinetic properties of the curing process can be determined (Figure 3, right). In this case, a first-order reaction with a time-constant of 370 ms was found. However, the exact nature of the reaction is likely dependent on the sample thickness and irradiation power, which are parameters that were not fully controlled for these measurements.

Faster measurements are certainly possible from the point of view of the spectrometer, however, the present setup appears to be limited by the lamp turn-on time and brightness. A different source (e.g. a laser) would allow for significantly faster measurements.

Conclusion

Rapid measurements with high signal to noise ratios were performed on a sample of glue during curing, giving insight on processes that happen on a millisecond timescale. This effectively bridges the previously-demonstrated capabilities of the IRis-F1 to measure on the microsecond timescale, with the second to minute timescale. It enables processes over many orders of magnitude in time to be studied. [1]


图3：(左) 以25ms的时间分辨率测量的固化过程的差异光谱；(右) 绝对谱图中两点的动力学轨迹
Figure 3: (left) Difference spectra of the curing process measured with a time-resolution of 25 ms; (right) kinetic traces at two points in the absolute spectrum